



STE13 Antibody, Biotin conjugated

Product Code	CSB-PA339395LD01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33894
Immunogen	Recombinant Saccharomyces cerevisiae Dipeptidyl aminopeptidase A protein (1-119AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidyl aminopeptidase A, DPAP A, YSCIV, STE13, YCI1, YOR219C, YOR50-9
Immunogen Species	Saccharomyces cerevisiae
Target Names	STE13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.