



KANK1 Antibody, FITC conjugated

Product Code	CSB-PA614803LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14678
Immunogen	Recombinant Human KN motif and ankyrin repeat domain-containing protein 1 protein (66-255AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KN motif and ankyrin repeat domain-containing protein 1, Ankyrin repeat domain-containing protein 15, Kidney ankyrin repeat-containing protein, KANK1, ANKRD15, KANK, KIAA0172
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	KANK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.