



ARHGEF6 Antibody, FITC conjugated

Product Code	CSB-PA618974LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15052
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 6 protein (108-239AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 6, Alpha-Pix, COOL-2, PAK-interacting exchange factor alpha, Rac/Cdc42 guanine nucleotide exchange factor 6, ARHGEF6, COOL2, KIAA0006, PIXA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	ARHGEF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.