



ARHGEF6 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA618974LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q15052 |
| Immunogen | Recombinant Human Rho guanine nucleotide exchange factor 6 protein (108-239AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Rho guanine nucleotide exchange factor 6, Alpha-Pix, COOL-2, PAK-interacting exchange factor alpha, Rac/Cdc42 guanine nucleotide exchange factor 6, ARHGEF6, COOL2, KIAA0006, PIXA |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal transduction |
| Target Names | ARHGEF6 |