



DUSP5 Antibody, FITC conjugated

Product Code	CSB-PA620976LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16690
Immunogen	Recombinant Human Dual specificity protein phosphatase 5 protein (56-181AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 5, Dual specificity protein phosphatase hVH3, DUSP5, VH3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	DUSP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.