



DUSP5 Antibody, Biotin conjugated

Product Code	CSB-PA620976LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16690
Immunogen	Recombinant Human Dual specificity protein phosphatase 5 protein (56-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 5, Dual specificity protein phosphatase hVH3, DUSP5, VH3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	DUSP5