




MSS51 Antibody

Product Code	CSB-PA681195LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4VC12
Immunogen	Recombinant Human Putative protein MSS51 homolog, mitochondrial protein (258-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative protein MSS51 homolog, mitochondrial, Zinc finger MYND domain-containing protein 17, MSS51, ZMYND17
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	MSS51

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, HepG2 whole cell lysate, A549 whole cell lysate
 All lanes: MSS51 antibody at 4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 52, 27, 11 kDa
 Observed band size: 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.