



MSS51 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA681195LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q4VC12 |
| Immunogen | Recombinant Human Putative protein MSS51 homolog, mitochondrial protein (258-365AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Putative protein MSS51 homolog, mitochondrial, Zinc finger MYND domain-containing protein 17, MSS51, ZMYND17 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell biology |
| Target Names | MSS51 |