



LEMD1 Antibody, HRP conjugated

Product Code	CSB-PA715035LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68G75
Immunogen	Recombinant Human LEM domain-containing protein 1 protein (44-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LEM domain-containing protein 1, Cancer/testis antigen 50, CT50, LEM domain protein 1, LEMP-1, LEMD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	LEMD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.