



cyp26b1 Antibody, Biotin conjugated

Product Code	CSB-PA721935LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6EIG3
Immunogen	Recombinant Zebrafish Cytochrome P450 26B1 protein (366-511AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 26B1, Retinoic acid-metabolizing cytochrome, cyp26b1
Immunogen Species	Zebrafish
Research Area	Cardiovascular
Target Names	cyp26b1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.