



LRIG3 Antibody, FITC conjugated

Product Code	CSB-PA751017LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXM1
Immunogen	Recombinant Human Leucine-rich repeats and immunoglobulin-like domains protein 3 protein (835-985AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeats and immunoglobulin-like domains protein 3, LIG-3, LRIG3, LIG3, UNQ287/PRO326/PRO335
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LRIG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.