



dye Antibody, Biotin conjugated

Product Code	CSB-PA772049LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7ZU29
Immunogen	Recombinant Zebrafish Nuclear pore complex protein Nup93 protein (1-109AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup93, 93 kDa nucleoporin, Dead eye protein, Nucleoporin Nup93, dye
Immunogen Species	Zebrafish
Target Names	dye