



dye Antibody, Biotin conjugated

Product Code	CSB-PA772049LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7ZU29
Immunogen	Recombinant Zebrafish Nuclear pore complex protein Nup93 protein (1-109AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup93, 93 kDa nucleoporin, Dead eye protein, Nucleoporin Nup93, dye
Immunogen Species	Zebrafish
Target Names	dye
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.