



CLEC14A Antibody, HRP conjugated

Product Code	CSB-PA801222LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86T13
Immunogen	Recombinant Human C-type lectin domain family 14 member A protein (278-349AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 14 member A, Epidermal growth factor receptor 5, EGFR-5, CLEC14A, C14orf27, EGFR5, UNQ236/PRO269
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	CLEC14A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.