



TIAM2 Antibody, HRP conjugated

Product Code	CSB-PA818235LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVF5
Immunogen	Recombinant Human T-lymphoma invasion and metastasis-inducing protein 2 protein (1500-1627AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-lymphoma invasion and metastasis-inducing protein 2, TIAM-2, SIF and TIAM1-like exchange factor, TIAM2, KIAA2016, STEF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	TIAM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.