



TIAM2 Antibody, Biotin conjugated

Product Code	CSB-PA818235LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVF5
Immunogen	Recombinant Human T-lymphoma invasion and metastasis-inducing protein 2 protein (1500-1627AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-lymphoma invasion and metastasis-inducing protein 2, TIAM-2, SIF and TIAM1-like exchange factor, TIAM2, KIAA2016, STEF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	TIAM2