



PIAS4 Antibody, HRP conjugated

Product Code	CSB-PA822692LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2W9
Immunogen	Recombinant Human E3 SUMO-protein ligase PIAS4 protein (66-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 SUMO-protein ligase PIAS4, PIASy, Protein inhibitor of activated STAT protein 4, Protein inhibitor of activated STAT protein gamma, PIAS-gamma, RING-type E3 ubiquitin transferase PIAS4, PIAS4, PIASG
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PIAS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.