



NSUN5 Antibody, Biotin conjugated

Product Code	CSB-PA839386LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P11
Immunogen	Recombinant Human Probable 28S rRNA (cytosine-C(5))-methyltransferase protein (362-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable 28S rRNA (cytosine-C(5))-methyltransferase, NOL1-related protein, NOL1R, NOL1/NOP2/Sun domain family member 5, Williams-Beuren syndrome chromosomal region 20A protein, NSUN5, NSUN5A, WBSCR20, WBSCR20A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	NSUN5