



NSUN5 Antibody, Biotin conjugated

Product Code	CSB-PA839386LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P11
Immunogen	Recombinant Human Probable 28S rRNA (cytosine-C(5))-methyltransferase protein (362-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable 28S rRNA (cytosine-C(5))-methyltransferase, NOL1-related protein, NOL1R, NOL1/NOP2/Sun domain family member 5, Williams-Beuren syndrome chromosomal region 20A protein, NSUN5, NSUN5A, WBSCR20, WBSCR20A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	NSUN5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.