



FIZ1 Antibody, Biotin conjugated

Product Code	CSB-PA839415LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SL8
Immunogen	Recombinant Human Flt3-interacting zinc finger protein 1 protein (252-410AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flt3-interacting zinc finger protein 1, Zinc finger protein 798, FIZ1, ZNF798
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FIZ1