




TSPEAR Antibody

Product Code	CSB-PA848817LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WU66
Immunogen	Recombinant Human Thrombospondin-type laminin G domain and EAR repeat-containing protein (248-433AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombospondin-type laminin G domain and EAR repeat-containing protein, TSP-EAR, TSPEAR, C21orf29
Immunogen Species	Homo sapiens (Human)
Target Names	TSPEAR

Image



Western Blot

Positive WB detected in: Raji whole cell lysate

All lanes: TSPEAR antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75, 67 kDa

Observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.