



SLCO2A1 Antibody, Biotin conjugated

Product Code	CSB-PA856440LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92959
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 2A1 protein (416-518AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier organic anion transporter family member 2A1, Prostaglandin transporter, PGT, Solute carrier family 21 member 2, SLCO2A1, OATP2A1, SLC21A2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	SLCO2A1