




GALNT9 Antibody

Product Code	CSB-PA862072LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCQ5
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 9 protein (184-338AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	GALNT9

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate, 293T whole cell lysate

All lanes: GALNT9 antibody at 2.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 69, 28 kDa

Observed band size: 69 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.