



GALNT9 Antibody, Biotin conjugated

Product Code	CSB-PA862072LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCQ5
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 9 protein (184-338AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polypeptide N-acetylgalactosaminyltransferase 9, Polypeptide GalNAc transferase 9, GalNAc-T9, pp-GaNTase 9, Protein-UDP acetylgalactosaminyltransferase 9, polypeptide N-acetylgalactosaminyltransferase 9, GALNT9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	GALNT9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.