



mycn Antibody, Biotin conjugated

Product Code	CSB-PA868688LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9PSJ0
Immunogen	Recombinant Zebrafish N-myc protein (1-93AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-myc protein, zN-Myc, mycn, nmyc1
Immunogen Species	Zebrafish
Target Names	mycn