



MRO Antibody, Biotin conjugated

Product Code	CSB-PA874842LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYG7
Immunogen	Recombinant Human Protein maestro protein (12-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein maestro, Male-specific transcription in the developing reproductive organs, Protein B29, MRO, B29, C18orf3
Immunogen Species	Homo sapiens (Human)
Target Names	MRO