



MRO Antibody, Biotin conjugated

Product Code	CSB-PA874842LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYG7
Immunogen	Recombinant Human Protein maestro protein (12-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein maestro, Male-specific transcription in the developing reproductive organs, Protein B29, MRO, B29, C18orf3
Immunogen Species	Homo sapiens (Human)
Target Names	MRO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.