



YY1AP1 Antibody, HRP conjugated

Product Code	CSB-PA884488LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H869
Immunogen	Recombinant Human YY1-associated protein 1 protein (335-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	YY1-associated protein 1, Hepatocellular carcinoma susceptibility protein1, YY1AP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	YY1AP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.