



PURG Antibody, Biotin conjugated

Product Code	CSB-PA891529LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJV8
Immunogen	Recombinant Human Purine-rich element-binding protein gamma protein (242-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Purine-rich element-binding protein gamma, PURG
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PURG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.