



RANBP3L Antibody, Biotin conjugated

Product Code	CSB-PA803130LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VV4
Immunogen	Recombinant Human Lck-interacting transmembrane adapter 1 protein (170-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ran-binding protein 3-like, RANBP3L
Immunogen Species	Homo sapiens (Human)
Target Names	RANBP3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.