




CCDC170 Antibody

Product Code	CSB-PA815585LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYT3
Immunogen	Recombinant Human Coiled-coil domain-containing protein 170 protein (574-651AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 170, CCDC170, C6orf97
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	CCDC170

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, PC-3 whole cell lysate

All lanes: CCDC170 antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 83 kDa

Observed band size: 83 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.