



RHBDF1 Antibody, HRP conjugated

Product Code	CSB-PA856909LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CC6
Immunogen	Recombinant Human Complement factor H-related protein 1 protein (132-221AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive rhomboid protein 1, iRhom1, Epidermal growth factor receptor-related protein, Rhomboid 5 homolog 1, Rhomboid family member 1, p100hRho, RHBDF1, C16orf8, DIST1, IRHOM1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	RHBDF1