



# RHBDF1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA856909LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96CC6
<b>Immunogen</b>	Recombinant Human Complement factor H-related protein 1 protein (132-221AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Inactive rhomboid protein 1, iRhom1, Epidermal growth factor receptor-related protein, Rhomboid 5 homolog 1, Rhomboid family member 1, p100hRho, RHBDF1, C16orf8, DIST1, IRHOM1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal transduction
<b>Target Names</b>	RHBDF1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.