



TIMMDC1 Antibody, HRP conjugated

Product Code	CSB-PA885685LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPL8
Immunogen	Recombinant Human E3 ubiquitin-protein ligase UHRF2 protein (444-585AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complex I assembly factor TIMMDC1, mitochondrial, Protein M5-14, Translocase of inner mitochondrial membrane domain-containing protein 1, TIMM domain containing-protein 1, TIMMDC1, C3orf1, UNQ247/PRO284
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	TIMMDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.