



CD47 Antibody, Biotin conjugated

Product Code	CSB-PA004940LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08722
Immunogen	Recombinant Human Leukocyte surface antigen CD47 protein (44-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte surface antigen CD47, Antigenic surface determinant protein OA3, Integrin-associated protein, IAP, Protein MER6, , CD47, CD47, MER6
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer; Immunology
Target Names	CD47
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.