



COL4A5 Antibody, HRP conjugated

Product Code	CSB-PA005746LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29400
Immunogen	Recombinant Human Collagen alpha-5(IV) chain protein (1184-1235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-5(IV) chain, COL4A5
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer; Signal transduction; Stem cells
Target Names	COL4A5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.