



ENTPD1 Antibody, Biotin conjugated

Product Code	CSB-PA007690LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49961
Immunogen	Recombinant Human Ectonucleoside triphosphate diphosphohydrolase 1 protein (316-474AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ectonucleoside triphosphate diphosphohydrolase 1, NTPDase 1, Ecto-ATP diphosphohydrolase 1, Ecto-ATPDase 1, Ecto-ATPase 1, Ecto-apyrase, Lymphoid cell activation antigen, , CD39, ENTPD1, CD39
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ENTPD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.