



GLI4 Antibody, HRP conjugated

Product Code	CSB-PA009502LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10075
Immunogen	Recombinant Human Zinc finger protein GLI4 protein (7-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein GLI4, Krueppel-related zinc finger protein 4, Protein HKR4, GLI4, HKR4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GLI4