



GLI4 Antibody, Biotin conjugated

Product Code	CSB-PA009502LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10075
Immunogen	Recombinant Human Zinc finger protein GLI4 protein (7-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein GLI4, Krueppel-related zinc finger protein 4, Protein HKR4, GLI4, HKR4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GLI4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.