



GRM6 Antibody, Biotin conjugated

Product Code	CSB-PA009936LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15303
Immunogen	Recombinant Human Metabotropic glutamate receptor 6 protein (370-531AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metabotropic glutamate receptor 6, mGluR6, GRM6, GPRC1F, MGLUR6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience; Signal transduction
Target Names	GRM6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.