



HOXC6 Antibody, Biotin conjugated

Product Code	CSB-PA010677LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09630
Immunogen	Recombinant Human Homeobox protein Hox-C6 protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-C6, Homeobox protein CP25, Homeobox protein HHO.C8, Homeobox protein Hox-3C, HOXC6, HOX3C
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.