



HRG Antibody, Biotin conjugated

Product Code	CSB-PA010736LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04196
Immunogen	Recombinant Human Histidine-rich glycoprotein protein (412-511AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histidine-rich glycoprotein, Histidine-proline-rich glycoprotein, HPRG, HRG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HRG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.