



IMPG2 Antibody, Biotin conjugated

Product Code	CSB-PA011704LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZV3
Immunogen	Recombinant Human Interphotoreceptor matrix proteoglycan 2 protein (616-759AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interphotoreceptor matrix proteoglycan 2, Interphotoreceptor matrix proteoglycan of 200 kDa, IPM 200, Sialoprotein associated with cones and rods proteoglycan, Spacrcan, IMPG2, IPM200
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	IMPG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.