



KRT16 Antibody, Biotin conjugated

Product Code	CSB-PA012515LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08779
Immunogen	Recombinant Human Keratin, type I cytoskeletal 16 protein (265-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 16, Cytokeratin-16, CK-16, Keratin-16, K16, KRT16, KRT16A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	KRT16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.