




MATN3 Antibody

Product Code	CSB-PA013522EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15232
Immunogen	Recombinant Human Matrilin-3 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrilin-3, MATN3
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental biology; Signal transduction
Target Names	MATN3

Image



Western Blot

Positive WB detected in: RAW264.7 whole cell lysate, SH-SY5Y whole cell lysate, HepG2 whole cell lysate, HeLa whole cell lysate, MCF-7 whole cell lysate

All lanes: MATN3 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53, 49 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.