



NOS2 Antibody, Biotin conjugated

Product Code	CSB-PA015943LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35228
Immunogen	Recombinant Human Nitric oxide synthase, inducible protein (17-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitric oxide synthase, inducible, Hepatocyte NOS, HEP-NOS, Inducible NO synthase, Inducible NOS, iNOS, NOS type II, Peptidyl-cysteine S-nitrosylase NOS2, NOS2, NOS2A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience; Cancer; Metabolism
Target Names	NOS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.