



NRL Antibody, HRP conjugated

Product Code	CSB-PA016086LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54845
Immunogen	Recombinant Human Neural retina-specific leucine zipper protein (24-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neural retina-specific leucine zipper protein, NRL, NRL, D14S46E
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling; Neuroscience
Target Names	NRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.