



NUP62 Antibody, FITC conjugated

| | |
|----------------------------|---------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA016204LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P37198 |
| Immunogen | Recombinant Human Nuclear pore glycoprotein p62 protein (173-522AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nuclear pore glycoprotein p62, 62 kDa nucleoporin, Nucleoporin Nup62, NUP62 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling; Cancer; Cell biology; Tags & Cell Markers; Signal transduction |
| Target Names | NUP62 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |