



PRPF40A Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA018769LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75400 |
| Immunogen | Recombinant Human Pre-mRNA-processing factor 40 homolog A protein (120-266AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pre-mRNA-processing factor 40 homolog A, Fas ligand-associated factor 1, Formin-binding protein 11, Formin-binding protein 3, Huntingtin yeast partner A, Huntingtin-interacting protein 10, HIP-10, Huntingtin-interacting protein A, Renal carcinoma antigen NY-REN-6, PRPF40A, FBP11, FLAF1, FNBP3, HIP10, HYPA, HSPC225 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | PRPF40A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |