



PRPF40A Antibody, Biotin conjugated

Product Code	CSB-PA018769LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75400
Immunogen	Recombinant Human Pre-mRNA-processing factor 40 homolog A protein (120-266AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA-processing factor 40 homolog A, Fas ligand-associated factor 1, Formin-binding protein 11, Formin-binding protein 3, Huntingtin yeast partner A, Huntingtin-interacting protein 10, HIP-10, Huntingtin-interacting protein A, Renal carcinoma antigen NY-REN-6, PRPF40A, FBP11, FLAF1, FNBP3, HIP10, HYPA, HSPC225
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRPF40A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.